

clinell®

www.clinell.com

EFFECTIVE
AGAINST
BACTERIA &
VIRUSES



CLINELL HAND SANITISING LIQUID

Unique liquid based formula

Specifically designed for use within the healthcare sector to kill germs and instantly sanitise hands.

Dermatologically tested

The hand sanitising liquids contain moisturiser to promote healthy, soft skin.

Available in 2 spray bottle sizes

The 60ml liquids small and convenient size makes it ideal for carrying with you at all times – at home, in the office, in your car, in your handbag or when travelling. Also available in a large 520ml pump dispenser.

Cleans and refreshes

With natural extracts of aloe vera.

HAND SANITISING LIQUIDS

Sanitising liquids with unique formula for instantly disinfecting hands without the need for soap and water. Containing a unique blend of 75% alcohol and moisturisers providing a superior level of effectiveness against germs without drying hands out.




HOW TO CLEAN HANDS EFFECTIVELY




1. Squirt into the palm of the hand.



2. Rub hands palm to palm.




3. Rub back of each hand with the palm of other hand with fingers interlaced.




4. Rub with backs of fingers to opposing palms with fingers interlocked.




5. Rub palm to palm with fingers interlaced.



6. Rub each thumb with the other hand, in a circle.



7. Rub tips of fingers in opposite palm in a circular motion.



8. Rub each wrist with opposite hand. Leave to air dry.

Clinell Hand Sanitising Liquids contain a unique blend of 75% alcohol and additional moisturisers, providing a superior level of effectiveness against germs. The liquid based formulation is easy to apply through the spray dispenser and helps to prevent hands from feeling sticky after application. Provides hand disinfection when soap and water is unavailable.

The formulation contains a unique blend of moisturisers which helps to condition the skin, preventing dry, cracked hands often associated with other alcohol hand sanitisers, leaving your hands feeling naturally refreshed.

BACTERIA	CONTACT TIME	TEST METHOD
<i>A. baumannii</i>	10 sec	EN 13727
<i>Enterococcus faecium</i> (VRE)	10 sec	EN 13727
<i>E.coli</i>	10 sec	EN 13727
<i>Klebsiella pneumoniae</i> (ESBL)	10 sec	EN 13727
<i>Pseudomonas aeruginosa</i>	10 sec	EN 13727
MRSA	10 sec	EN 13727
<i>Staphylococcus aureus</i>	10 sec	EN 13727
<i>Streptococcus pneumoniae</i>	10 sec	EN 13727
<i>K.pneumoniae</i> (CPE)	10 sec	EN 13727
<i>Enterobacter spp.</i>	10 sec	EN 13727
<i>Enterococcus hirae</i>	10 sec	EN 13727
<i>Carbapenem resistant enterobacteriaceae</i> (CRE)	10 sec	EN 13727
YEAST/FUNGI		
<i>C. albicans</i>	10 sec	EN 13624
<i>C. auris</i>	10 sec	EN 13624
<i>A. brasiliensis</i>	30 sec	EN 13624
MYCOBACTERIA		
<i>M. terrae</i>	15 sec	EN 14348
<i>M. avium</i>	15 sec	EN 14348
<i>M. smegmatis</i>	15 sec	EN 14348
VIRUSES		
Norovirus	15 sec	EN 14476
Influenza	15 sec	EN 14476
HIV	15 sec	EN 14476
Hepatitis C	15 sec	EN 14476
Vaccinia virus	30 sec	EN 14476
Adenovirus	30 sec	EN 14476
Murine coronavirus (SARS-CoV-2 surrogate)	30 sec	EN 14476
HAND ANTISEPTIC TESTING		
Hygienic handrub	30 sec	EN 1500
Surgical hand disinfection	90 sec	EN 12791

PRODUCT	UNIT OF ISSUE	UNITS PER BOX	CODE
Hand Sanitising Liquid	1 x 60ml	96	CFA60UK
Hand Sanitising Liquid pump with head	1 x 520ml	10	CFA520UK

Use biocides safely. Always read the label and product information before use.

GAMA Healthcare Ltd., 2 Regal Way, Watford, Hertfordshire, WD24 4YJ.
T: +44 (0)20 7993 0030 E: info@gamahealthcare.com www.gamahealthcare.com

gama
healthcare